## Quick Guide

# ImmuView<sup>®</sup> Reader 2.0

English (UK)

## Introduction

This Quick Guide provides an overview of how to operate the ImmuView® Reader 2.0 with SSI Diagnostica AVS approved immunochromatographic *in vitro* diagnostics tests.

For comprehensive instructions refer to the ImmuView® Reader 2.0 User Manual available online at: <https://ssidiagnostica.com/ifu/immuview-reader/> or by scanning the QR code

## Intended purpose

The ImmuView® Reader 2.0 is an adjunctive instrument for automated reading and interpretation of results for immunochromatographic *in vitro* diagnostics tests.

The instrument is to be used by laboratory professionals or health care professionals only.

## Instrument overview

Figure 1: Front view


Figure 2: Rear view

## Starting the instrument

To power on the instrument, press and hold the power button (located on the left side, Figure 2) for 2–3 seconds until an audible tone is heard.

The startup process includes an automated self-test to verify internal system functionality and takes a few minutes to complete.

## Navigation bar icons

**Test** – Initiate a new test

**Export** – Export results (available if connected to LIS or if a FAT32 formatted USB key is present in one of the USB ports)

**Home** – Return to the main menu

**Print** – Print results (requires connected printer)

## Home screen (main menu)


*NOTE. Screen images shown are for illustrative purposes only. Test availability varies depending on the instrument configuration and availability in the market.*

Upon completion of the startup sequence, the Home menu appears.

Enter your login credentials when prompted for a username and password.


Available Options:

- TEST – Initiate a patient sample, QC test, or instrument check test
- RESULTS – View previously recorded test results
- SETTINGS – Adjust instrument configuration

## Performing an Instrument Check test

- From the Home menu, select "TEST", then choose "Instrument Check."
- Insert the Instrument Check Cartridge into the designated recess, ensuring the arrow faces upward and toward the instrument.


- The test begins automatically upon closing the Cartridge Drawer.
  - A "Pass" result means that the Instrument Check has passed.
  - A "Fail" result means that the Instrument Check has failed.
  - An "Overdue" result means that the Instrument Check pass status validity period has ended, or an Instrument Check has not previously been run.

*NOTE. The instrument is configured to prompt for an Instrument Check test every 7 days. If not completed or if the test fails, sample testing will be locked out until a valid Instrument Check test is performed.*

## Running a patient sample test


Consult the test-specific IFU and Quick Guide for instructions on how to perform the test.

### Step 1: Select "TEST" from the Home menu


### Step 2: Choose test mode

- Select "Perform Test" for immediate reading.
- Select "Incubate & Read" for timed incubation and reading.


*Note: "Incubate & Read" may not be available on all configurations or for all test types.*

### Step 3: Choose type of test


- If required, input or scan the test kit lot information.
- Confirm by pressing ✓ (Accept).


### Step 4: Enter patient ID

- Input or scan a unique patient ID (avoid using patient names).
- Confirm by pressing ✓ (Accept).

### Step 5: Insert test

Open the Cartridge Drawer and place the test in the designated recess (refer to the test-specific Quick Guide for details).



*NOTE. Be attentive to the orientation of the test when placing in the Cartridge Drawer - incorrect orientation may result in an invalid result.*

## Step 6:

### Close Cartridge Drawer


- In “Immediate Read” mode (Perform Test), the test is analysed immediately.
- In “Incubate & Read” mode, a countdown timer begins, and the test is analysed throughout the incubation period.

*NOTE. Ensure the Cartridge Drawer remains closed during analysis.*

### Viewing results

- Upon completion, an audible tone signals that results are ready.
- Results are displayed and saved automatically.
- If the result is “Invalid,” refer to the assay-specific IFU or Quick Guide.
- You may now safely remove the test.
- To review saved results, select “RESULTS” from the Home menu.

## Quality certification

SSI Diagnostica A/S ensures that the development, production, and distribution of *in vitro* diagnostics comply with EN ISO 13485 standard.

REF  
213414

IVD

## Contact & ordering information

### Legal Manufacturer

SSI Diagnostica A/S  
Herredsvejen 2  
3400 Hillerød  
Denmark  
T +45 4829 9100  
info@ssidiagnostica.com  
ssidiagnostica.com

### UK Representative (UKRP)

Emergo Consulting (UK) Limited  
c/o Cr360 – UL International  
Compass House, Vision Park Histon  
Cambridge CB24 9BZ  
United Kingdom

### Swiss Representative

CH REP

MedEnvoy Switzerland  
Gotthardstrasse 28  
6302 Zug  
Switzerland

© 2025 SSI Diagnostica Group. All rights reserved. ImmuView and the ImmuView logo are trademarks of SSI Diagnostica A/S. C. DIFF QUIK CHEK COMPLETE is a trademark of TECHLAB, Inc., under license.

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