Quick Guide

How to use the
ImmuView[®] Reader 2.0 with the
C. DIFF QUIK CHEK COMPLETE[®] test

English (UK)

Introduction

This Quick Guide provides an overview of how to operate the ImmuView® Reader 2.0 with the TECHLAB® *C. DIFF QUIK CHEK COMPLETE*® test.

For additional instructions, refer to:


- ImmuView® Reader 2.0 User Manual available at <https://ssidiagnostica.com/ifu/immuview-reader/> or by scanning the QR code
- ImmuView® Reader 2.0 Quick Guide
- *C. DIFF QUIK CHEK COMPLETE*® Instructions for Use

Running a patient sample test


Perform the test as described in the *C. DIFF QUIK CHEK COMPLETE®* test Instructions for Use.

NOTE. Screen images shown are for illustrative purposes only. Test availability varies depending on the instrument configuration and availability in the market.

Step 1:
Select "TEST" from the Home menu

Step 2:
Select "Perform Test"


Step 3:

Select "C. DIFF QUIK CHEK COMPLETE"


- If required, input or scan the test kit lot information.
- Confirm by pressing ✓ (Accept).

Step 4:

Enter patient ID

- Input or scan a unique patient ID (avoid using patient names).
- Confirm by pressing ✓ (Accept).


Step 5:

Place Membrane Device in Cartridge Drawer

Open the Cartridge Drawer and place the Membrane Device into the indicated recess as shown, with the Sample Well in the upper-right corner.


NOTE. Be attentive to the orientation of the test when placing in the Cartridge Drawer - incorrect orientation may result in an invalid result.

Step 6:
Close the Cartridge Drawer

The test is analysed immediately.

NOTE. Ensure the Cartridge Drawer remains closed during analysis.

Viewing results

- Upon completion, an audible tone signals that results are ready.
- Results are displayed and saved automatically.
- If the result is "Invalid," refer to the IFU and retest the sample.
- You may now safely remove the Membrane Device.

Quality certification

SSI Diagnostica A/S ensures that the development, production, and distribution of *in vitro* diagnostics comply with EN ISO 13485 standard.

REF
213414

IVD

ImmuView® Reader 2.0 contact & ordering information

Legal Manufacturer

SSI Diagnostica A/S
Herredsvæjen 2
3400 Hillerød
Denmark
T +45 4829 9100
info@ssidiagnostica.com
ssidiagnostica.com

UK Representative (UKRP)

Emergo Consulting (UK) Limited
c/o Cr360 – UL International
Compass House, Vision Park Histon
Cambridge CB24 9BZ
United Kingdom

Swiss Representative

CH REP

MedEnvoy Switzerland
Gotthardstrasse 28
6302 Zug
Switzerland

© 2025 SSI Diagnostica Group. All rights reserved. ImmuView and the ImmuView logo are trademarks of SSI Diagnostica A/S. C. DIFF QUIK CHEK COMPLETE is a trademark of TECHLAB, Inc., under license.

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