

DIAGNOSE WITH CONFIDENCE

IMMUVIEW[®]

S. pneumoniae
and *Legionella*

Identify both with
just one test

For urine and CSF samples




SSI Diagnostica A/S
Herredsvej 2
3400 Hillerød
P +45 4829 9100

IMMUVIEW® DETECT MORE THAN SEROGROUP 1 OF *LEGIONELLA*

ImmuView® *Legionella* serogroups - control antigens¹

Positive	1 (All monoclonal subtypes)
	2, 3, 4, 5, 6, 7, 11, 12, 13, 14, 15
Negative	8, 9, 10



Legionella LOD for different rapid antigen tests²

Serogroup (subtypes)	ImmuView®	BinaxNOW™
SG1 (Philadelphia 1st)	10 ⁻⁵ x 1:8	10 ⁻⁵ x 1:8
SG1 (OLDA/Oxford)	10 ⁻⁵ x 1:4	10 ⁻⁵ x 1:4
SG1 (Allentown/France)	10 ⁻⁵ x 1:4	10 ⁻⁵ x 1:4
SG1 (Benidorm)	10 ⁻⁴ x 1:8	10 ⁻⁴ x 1:16 ^A
SG1 (Bellingham)	10 ⁻⁴ x 1:8	10 ⁻⁴ x 1:4 ^A
SG3	UD x 1:32	Negative^B
SG6	UD x 1:64	Negative^B
SG8	UD x 1:256	Negative^B
SG10	UD x 1:32	Negative^B

UD = Undiluted samples; A) Very faint lines; B) Negative when tested in undiluted samples

1. Lindsay et. Al, Evaluation of the rapid lateral flow ImmuView® urinary antigen test for the detection of *L. pneumophila* and *Streptococcus pneumoniae*, poster presented at ESGLI 2017.

2. Data adapted from; Andersen, et. Al, Sensitivity for *S. pneumoniae* and *L. pneumophila* in Different Rapid Urinary Antigen Assays; poster at ASM 2019

Product details
ImmuView® *S. pneumoniae* and *L. pneumophila* Urinary Antigen Test, Article No. 95389 . Kit Size: 22 tests.



DIAGNOSE WITH CONFIDENCE